

Product Datasheet

Mouse IL31 ELISA Kit (orb408948)

Description	Mouse Interleukin 31(IL-31) ELISA Kit
Reactivity	Mouse
Range	31.25 pg/mL-2000 pg/mL
Target	IL31
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, cell culture supernates, tissue homogenates
Uniprot ID	Q6EAL8
Sensitivity	7.81 pg/mL
Expiration Date	Please enquire.