

Product Datasheet

Mouse IL1F6 ELISA Kit (orb408952)

Description	Mouse Interleukin-36 alpha(IL1F6) ELISA kit
Reactivity	Mouse
Range	15.6 pg/mL-1000 pg/mL
Target	Il36a
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	Q9JLA2
Sensitivity	3.9 pg/mL
Expiration Date	Please enquire.