

Product Datasheet

Mouse IL1F8 ELISA Kit (orb408953)

Description	Mouse Interleukin-36 beta(IL1F8) ELISA kit
Reactivity	Mouse
Range	6.25 pg/mL-400 pg/mL
Target	IL1F8
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	Q9D6Z6
Sensitivity	1.56 pg/mL
Expiration Date	Please enquire.