

Product Datasheet

Mouse SOD3 ELISA Kit (orb408966)

Description	Mouse Extracellular superoxide dismutase [Cu-Zn](SOD3) ELISA kit
Reactivity	Mouse
Range	7.8 pg/mL-500 pg/mL
Target	SOD3
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	009164
Sensitivity	1.95 pg/mL
Expiration Date	Please enquire.