

Product Datasheet

Mouse ALT/GPT ELISA Kit (orb408971)

Description	Mouse Alanine Aminotransferase(ALT/GPT)ELISA Kit
Reactivity	Mouse
Range	0.156 mU/mL-10 mU/mL
Target	ALT/GPT
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, cell culture supernates, cell lysates
Uniprot ID	Q8QZR5
Sensitivity	0.039 mU/mL
Expiration Date	Please enquire.