

## Product Datasheet

### Mouse Complement 3 ELISA Kit (orb408975)

|                          |                                      |
|--------------------------|--------------------------------------|
| <b>Description</b>       | Mouse Complement 3, C3<br>ELISA Kit  |
| <b>Reactivity</b>        | Mouse                                |
| <b>Range</b>             | 93.75 ng/mL-6000 ng/mL               |
| <b>Target</b>            | C3                                   |
| <b>Note</b>              | For research use only                |
| <b>Application notes</b> | Detection Wavelength: 450<br>nm      |
| <b>Sample Types</b>      | serum, plasma, tissue<br>homogenates |
| <b>Uniprot ID</b>        | <b>P01027</b>                        |
| <b>Sensitivity</b>       | 23.44 ng/mL                          |
| <b>Expiration Date</b>   | Please enquire.                      |