

Product Datasheet

Mouse iFABP ELISA Kit (orb408986)

Description	Mouse intestinal fatty acid binding protein,iFABP ELISA Kit
Reactivity	Mouse
Range	6.25 pg/mL-400 pg/mL
Target	iFABP
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, cell culture supernates, tissue homogenates
Uniprot ID	P55050
Sensitivity	1.56 pg/mL
Expiration Date	Please enquire.