

## Product Datasheet

### Mouse TRH ELISA Kit (orb409001)

|                          |                                                   |
|--------------------------|---------------------------------------------------|
| <b>Description</b>       | Mouse thyrotropin-releasing hormone,TRH ELISA KIT |
| <b>Reactivity</b>        | Mouse                                             |
| <b>Range</b>             | 2.5 pg/mL-160 pg/mL                               |
| <b>Target</b>            | TRH                                               |
| <b>Note</b>              | For research use only                             |
| <b>Application notes</b> | Detection Wavelength: 450 nm                      |
| <b>Sample Types</b>      | serum                                             |
| <b>Uniprot ID</b>        | <b>Q62361</b>                                     |
| <b>Sensitivity</b>       | 0.625 pg/mL                                       |
| <b>Expiration Date</b>   | Please enquire.                                   |