

## Product Datasheet

### mouse FSH ELISA Kit (orb409002)

|                          |                                                  |
|--------------------------|--------------------------------------------------|
| <b>Description</b>       | mouse follicle-stimulating hormone,FSH ELISA Kit |
| <b>Reactivity</b>        | Mouse                                            |
| <b>Range</b>             | 4 mIU/ml-140 mIU/ml                              |
| <b>Target</b>            | FSH                                              |
| <b>Note</b>              | For research use only                            |
| <b>Application notes</b> | Detection Wavelength: 450 nm                     |
| <b>Sample Types</b>      | serum, plasma                                    |
| <b>Uniprot ID</b>        | <b>Q60687</b>                                    |
| <b>Sensitivity</b>       | 2.5 mIU/ml                                       |
| <b>Expiration Date</b>   | Please enquire.                                  |