

## Product Datasheet

### Mouse CRH ELISA Kit (orb409004)

|                          |                                                     |
|--------------------------|-----------------------------------------------------|
| <b>Description</b>       | Mouse Corticotropin-releasing Factor(CRF) ELISA Kit |
| <b>Reactivity</b>        | Mouse                                               |
| <b>Range</b>             | 2.5 pg/mL-40 pg/mL                                  |
| <b>Target</b>            | CRH                                                 |
| <b>Note</b>              | For research use only                               |
| <b>Application notes</b> | Detection Wavelength: 450 nm                        |
| <b>Sample Types</b>      | serum, plasma                                       |
| <b>Uniprot ID</b>        | <b>Q8CITO</b>                                       |
| <b>Sensitivity</b>       | 0.675 pg/mL                                         |
| <b>Expiration Date</b>   | Please enquire.                                     |