

Product Datasheet

Mouse MCP1 ELISA Kit (orb409010)

Description	Mouse monocyte chemotactic protein 1/monocyte chemotactic and activating factor,MCP-1/MCAF ELISA kit
Reactivity	Mouse
Range	78 pg/mL-5000 pg/mL
Target	CCL2
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	P10148
Sensitivity	19.5 pg/mL
Expiration Date	Please enquire.