

Product Datasheet

Mouse Chat ELISA Kit (orb409014)

Description	Mouse choline acetyltransferase(ChAT)ELISA Kit
Reactivity	Mouse
Range	15.6 pg/mL-1000 pg/mL
Target	CHAT
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates, cell culture supernates
Uniprot ID	Q03059
Sensitivity	3.9 pg/mL
Expiration Date	Please enquire.