

Product Datasheet

Mouse MCT ELISA Kit (orb409044)

Description	Mouse Mast Cell Tryptase(MCT) ELISA Kit
Reactivity	Mouse
Range	31.25 pg/mL-2000 pg/mL
Target	MCT
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, cell culture supernates, tissue homogenates
Uniprot ID	Q02844
Sensitivity	7.8 pg/mL
Expiration Date	Please enquire.