

Product Datasheet

Mouse HGFAC ELISA Kit (orb409054)

Description	Mouse Hepatocyte growth factor activator(HGFAC) ELISA kit
Reactivity	Mouse
Range	31.25 pg/mL-2000 pg/mL
Target	HGFAC
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, cell culture supernates, tissue homogenates
Uniprot ID	Q9R098
Sensitivity	7.8 pg/mL
Expiration Date	Please enquire.