

Product Datasheet

Mouse MET ELISA Kit (orb409055)

Description	Mouse Hepatocyte Growth Factor Receptor(C-MET/HGFR)ELISA Kit
Reactivity	Mouse
Range	0.078 ng/mL-5 ng/mL
Target	MET
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	P16056
Sensitivity	0.02 ng/mL
Expiration Date	Please enquire.