

Product Datasheet

Mouse IFNA2 ELISA Kit (orb409059)

Description	Mouse interferon alpha-2(IFNA2) ELISA kit
Reactivity	Mouse
Range	9.38 pg/mL-600 pg/mL
Target	IFNA2
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates, cell lysates
Uniprot ID	P01573
Sensitivity	2.35 pg/mL
Expiration Date	Please enquire.