

Product Datasheet

Mouse KITLG ELISA Kit (orb409063)

Description	Mouse Stem cell factor/mast cell growth factor, SCF/MGF ELISA kit
Reactivity	Mouse
Range	3.12 pg/mL-200 pg/mL
Target	KITLG
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	P20826
Sensitivity	0.78 pg/mL
Expiration Date	Please enquire.