

Product Datasheet

Mouse HMGB1 ELISA Kit (orb409067)

Description	Mouse High mobility group protein B1,HMGB-1 ELISA Kit
Reactivity	Mouse
Range	15.6 pg/mL-1000 pg/mL
Target	HMGB1
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates, cell lysates
Uniprot ID	P63158
Sensitivity	3.9 pg/mL
Expiration Date	Please enquire.