

Product Datasheet

Mouse GRIN2B ELISA Kit (orb409070)

Description	Mouse Glutamate [NMDA] receptor subunit epsilon-2 (GRIN2B) ELISA kit
Reactivity	Mouse
Range	31.25 pg/mL-2000 pg/mL
Target	GRIN2B
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates, cell lysates
Uniprot ID	Q01097
Sensitivity	7.8 pg/mL
Expiration Date	Please enquire.