

Product Datasheet

Mouse BMP2 ELISA Kit (orb409076)

Description	Mouse Bone morphogenetic protein 2,BMP-2 ELISA Kit
Reactivity	Mouse
Range	1.56 pg/mL-100 pg/mL
Target	BMP2
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, cell culture supernates, tissue homogenates
Uniprot ID	P21274
Sensitivity	0.39 pg/mL
Expiration Date	Please enquire.