

Product Datasheet

Mouse BMP7 ELISA Kit (orb409078)

Description	Mouse Bone morphogenetic protein 7,BMP-7 ELISA Kit
Reactivity	Mouse
Range	0.625 pg/mL-40 pg/mL
Target	BMP7
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, cell culture supernates, tissue homogenates
Uniprot ID	P23359
Sensitivity	0.156 pg/mL
Expiration Date	Please enquire.