

Product Datasheet

Mouse GDF2 ELISA Kit (orb409079)

Description	Mouse bone morphogenetic protein 9,BMP-9 ELISA Kit
Reactivity	Mouse
Range	0.312 pg/mL-20 pg/mL
Target	GDF2
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, cell culture supernates, tissue homogenates
Uniprot ID	Q9WV56
Sensitivity	0.078 pg/mL
Expiration Date	Please enquire.