

Product Datasheet

Mouse BGLAP ELISA Kit (orb409082)

Description	Mouse Osteocalcin/Bone gla protein,OT/BGP ELISA kit
Reactivity	Mouse
Range	78 pg/ml - 5000 pg/ml
Target	BGLAP
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	P86546
Sensitivity	19.5 pg/ml
Expiration Date	Please enquire.