

Product Datasheet

Mouse BMP6 ELISA Kit (orb409087)

Description	Mouse Bone morphogenetic protein 6,BMP-6 ELISA Kit
Reactivity	Mouse
Range	31.25 pg/mL-2000 pg/mL
Target	BMP6
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, cell culture supernates, tissue homogenates
Uniprot ID	P20722
Sensitivity	7.81 pg/mL
Expiration Date	Please enquire.