

Product Datasheet

Mouse PPARGC1A ELISA Kit (orb409096)

Description	Mouse Peroxisome proliferator-activated receptor gamma coactivator 1-alpha (PPARGC1A) ELISA kit
Reactivity	Mouse
Range	23.5 pg/mL-1500 pg/mL
Target	PPARGC1A
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	070343
Sensitivity	5.86 pg/mL
Expiration Date	Please enquire.