

Product Datasheet

Mouse PPARG ELISA Kit (orb409099)

Description	Mouse Peroxisome Proliferator-activated receptor gamma;;PPAR-gamma; ELISA Kit
Reactivity	Mouse
Range	0.625 pg/mL-40 pg/mL
Target	PPARG
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	P37238
Sensitivity	0.156 pg/mL
Expiration Date	Please enquire.