

Product Datasheet

Mouse RNASE2 ELISA Kit (orb409106)

Description	Mouse ribonuclease, RNase A family, 2 (liver, eosinophil-derived neurotoxin) (RNASE2) ELISA kit
Reactivity	Mouse
Range	3.12 ng/mL-200 ng/mL
Target	RNASE2
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	O35291
Sensitivity	0.78 ng/mL
Expiration Date	Please enquire.