

Product Datasheet

Mouse NFkB p65 ELISA Kit (orb409109)

Description	Mouse nuclear factor-kappaB p65,NF-kappaB p65 ELISA Kit
Reactivity	Mouse
Range	31.25 pg/mL-2000 pg/mL
Target	RELA
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	Q04207
Sensitivity	7.8 pg/mL
Expiration Date	Please enquire.