

Product Datasheet

Mouse Cox2 ELISA Kit (orb409114)

Description	Mouse cyclooxygenase-2, COX-2 ELISA Kit
Reactivity	Mouse
Range	31.25 pg/mL-2000 pg/mL
Target	PTGS2
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates, cell lysates
Uniprot ID	Q05769
Sensitivity	7.8 pg/mL
Expiration Date	Please enquire.