

Product Datasheet

Mouse XDH ELISA Kit (orb409115)

Description	Mouse Xanthine dehydrogenase/oxidase [Includes: Xanthine dehydrogenase(XDH) ELISA kit
Reactivity	Mouse
Range	7.8 pg/mL-500 pg/mL
Target	XDH
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	Q00519
Sensitivity	1.95 pg/mL
Expiration Date	Please enquire.