

## Product Datasheet

### Mouse CK-MB ELISA Kit (orb409125)

|                          |                                                            |
|--------------------------|------------------------------------------------------------|
| <b>Description</b>       | Mouse Creatine Kinase MB isoenzyme,CK-MB ELISA Kit         |
| <b>Reactivity</b>        | Mouse                                                      |
| <b>Range</b>             | 1.25 mU/mL-80 mU/mL                                        |
| <b>Target</b>            | CK-MB                                                      |
| <b>Note</b>              | For research use only                                      |
| <b>Application notes</b> | Detection Wavelength: 450 nm                               |
| <b>Sample Types</b>      | serum, plasma, cell culture supernates, tissue homogenates |
| <b>Sensitivity</b>       | 0.312 mU/mL                                                |
| <b>Expiration Date</b>   | Please enquire.                                            |