

Product Datasheet

Mouse MMP2 ELISA Kit (orb409133)

Description	Mouse matrix metalloproteinase 2/Gelatinase A,MMP-2 ELISA kit
Reactivity	Mouse
Range	15.6 pg/mL-1000 pg/mL
Target	MMP2
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, cell culture supernates, tissue homogenates
Uniprot ID	P33434
Sensitivity	3.9 pg/mL
Expiration Date	Please enquire.