

Product Datasheet

Mouse MMP7 ELISA Kit (orb409136)

Description	Mouse matrix metallopeptidase 7 (matrilysin, uterine) (MMP7) ELISA kit
Reactivity	Mouse
Range	31.25 pg/mL-2000 pg/mL
Target	MMP7
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, cell culture supernates, tissue homogenates
Uniprot ID	Q10738
Sensitivity	7.81 pg/mL
Expiration Date	Please enquire.