

Product Datasheet

Mouse MMP8 ELISA Kit (orb409137)

Description	Mouse matrix metalloproteinase 8/Neutrophil collagenase,MMP-8 ELISA kit
Reactivity	Mouse
Range	15.6 pg/mL-1000 pg/mL
Target	MMP8
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, cell culture supernates, tissue homogenates
Uniprot ID	O70138
Sensitivity	3.9 ng/mL
Expiration Date	Please enquire.