

## Product Datasheet

### Mouse TIMP2 ELISA Kit (orb409138)

|                          |                                                                 |
|--------------------------|-----------------------------------------------------------------|
| <b>Description</b>       | Mouse tissue inhibitors of metalloproteinase 2,TIMP-2 ELISA kit |
| <b>Reactivity</b>        | Mouse                                                           |
| <b>Range</b>             | 125 pg/mL-8000 pg/mL                                            |
| <b>Target</b>            | TIMP2                                                           |
| <b>Note</b>              | For research use only                                           |
| <b>Application notes</b> | Detection Wavelength: 450 nm                                    |
| <b>Sample Types</b>      | serum, plasma, cell culture supernates, tissue homogenates      |
| <b>Uniprot ID</b>        | <b>P25785</b>                                                   |
| <b>Sensitivity</b>       | 31.25 pg/mL                                                     |
| <b>Expiration Date</b>   | Please enquire.                                                 |