

Product Datasheet

Mouse TIMP3 ELISA Kit (orb409139)

Description	Mouse tissue inhibitors of metalloproteinase 3,TIMP-3 ELISA kit
Reactivity	Mouse
Range	46.88 pg/mL-3000 pg/mL
Target	TIMP3
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	P39876
Sensitivity	11.72 pg/mL
Expiration Date	Please enquire.