

Product Datasheet

Mouse SDF1 ELISA Kit (orb409142)

Description	Mouse Stromal cell derived factor 1 β ,SDF-1 β ELISA kit
Reactivity	Mouse
Range	31.25 pg/mL-2000 pg/mL
Target	CXCL12
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, cell culture supernates, tissue homogenates
Uniprot ID	P40224
Sensitivity	7.81 pg/mL
Expiration Date	Please enquire.