

Product Datasheet

Mouse COL18A1/ES ELISA Kit (orb409159)

Description	Mouse Collagen alpha-1(XVIII) chain (COL18A1/ES) ELISA kit
Reactivity	Mouse
Range	0.156 ng/mL-10 ng/mL
Target	COL18A1/ES
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, cell culture supernates, tissue homogenates
Uniprot ID	P39061
Sensitivity	0.039 ng/mL
Expiration Date	Please enquire.