

Product Datasheet

Mouse GDNF ELISA Kit (orb409160)

Description	Mouse glial cell line-derived neurotrophic factor,GDNF ELISA KIT
Reactivity	Mouse
Range	10.9 pg/ml - 8000 pg/ml
Target	GDNF
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	P48540
Sensitivity	10.9 pg/ml
Expiration Date	Please enquire.