

Product Datasheet

Mouse CNTF ELISA Kit (orb409167)

Description	Mouse Ciliary Neurotrophic Factor,CNTF ELISA Kit
Reactivity	Mouse
Range	6.25 pg/mL-400 pg/mL
Target	CNTF
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates, cell culture supernates
Uniprot ID	P51642
Sensitivity	1.56 pg/mL
Expiration Date	Please enquire.