

Product Datasheet

Mouse ADAM8 ELISA Kit (orb409168)

Description	Mouse A Disintegrin And Metalloprotease 8,ADAM8 ELISA Kit
Reactivity	Mouse
Range	62.5 pg/mL-4000 pg/mL
Target	ADAM8
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	Q05910
Sensitivity	15.6 pg/mL
Expiration Date	Please enquire.