

Product Datasheet

Mouse CD206 ELISA Kit (orb409173)

Description	Mouse Macrophage mannose receptor 1(MRC1) ELISA kit
Reactivity	Mouse
Range	0.625 ng/mL-40 ng/mL
Target	MRC1
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates, cell lysates
Uniprot ID	Q61830
Sensitivity	0.156 ng/mL
Expiration Date	Please enquire.