

Product Datasheet

Mouse CSF1 ELISA Kit (orb409174)

Description	Mouse Macrophage Colony-Stimulating Factor,M-CSF ELISA kit
Reactivity	Mouse
Range	7.81 pg/mL-500 pg/mL
Target	CSF1
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	P07141
Sensitivity	1.95 pg/mL
Expiration Date	Please enquire.