

Product Datasheet

Mouse CCL4 ELISA Kit (orb409176)

Description	Mouse Macrophage Inflammatory Protein 1 β ,MIP-1 β ELISA kit
Reactivity	Mouse
Range	12.5 pg/mL-800 pg/mL
Target	CCL4
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	P14097
Sensitivity	3.12 pg/mL
Expiration Date	Please enquire.