

Product Datasheet

Mouse PIGR ELISA Kit (orb409182)

Description	Mouse Polymeric immunoglobulin receptor(PIGR) ELISA kit
Reactivity	Mouse
Range	31.25 pg/mL-2000 pg/mL
Target	PIGR
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, cell culture supernates, tissue homogenates
Uniprot ID	070570
Sensitivity	7.81 pg/mL
Expiration Date	Please enquire.