

Product Datasheet

Mouse ACP5 ELISA Kit (orb409190)

Description	Mouse tartrate-resistant acid phosphatase 5b(TRACP-5b) ELISA Kit
Reactivity	Mouse
Range	0.078 mIU/mL-5 mIU/mL
Target	ACP5
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	Q05117
Sensitivity	0.02 mIU/mL
Expiration Date	Please enquire.