

Product Datasheet

Mouse RANKL ELISA Kit (orb409214)

Description	Mouse soluble receptor activator of nuclear factor-kB ligand,sRANKL ELISA Kit
Reactivity	Mouse
Range	12.5 pg/ml-800 pg/ml.
Target	TNFSF11
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates.
Uniprot ID	035235
Sensitivity	3.12 pg/ml.
Expiration Date	Please enquire.