

Product Datasheet

Mouse MYH1 ELISA Kit (orb409215)

Description	Mouse soluble myosin heavy chain 1,sMHC-1 ELISA Kit
Reactivity	Mouse
Range	15.6 pg/mL-1000 pg/mL
Target	MYH1
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, cell culture supernates, tissue homogenates
Uniprot ID	Q5SX40
Sensitivity	3.9 pg/mL
Expiration Date	Please enquire.