

Product Datasheet

Mouse LOX-1 ELISA Kit (orb409216)

Description	Mouse soluble Lectin-like Oxidized Low Density Lipoprotein Receptor-1(sLOX-1)ELISA Kit
Reactivity	Mouse
Range	78 pg/mL-5000 pg/mL
Target	OLR1
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	Q9EQ09
Sensitivity	19.5 pg/mL
Expiration Date	Please enquire.