

Product Datasheet

Mouse TREM1 ELISA Kit (orb409217)

Description	Mouse soluble Triggering Receptor Expresses on Myeloid Cells-1,sTREM-1 ELISA Kit
Reactivity	Mouse
Range	15.6 pg/mL-1000 pg/mL
Target	TREM1
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	Q9JKE2
Sensitivity	3.9 pg/mL
Expiration Date	Please enquire.